EN

Konelab™ / T Series CREATININE (Enzymatic)

REF 981845 4 x 60 ml

THIS PACKAGE INSERT IS APPLICABLE FOR USE OUTSIDE THE US. ANY REFERENCE TO THE KONELAB SYSTEMS ALSO REFERS TO THE T SERIES.

INTENDED USE

For the *in vitro* quantitative determination of creatinine concentration in human serum, plasma or urine on Konelab analyzers using enzymatic method. All test results must be interpreted with regard to the clinical context.


SUMMARY (1, 2)

Creatine is synthesized in the kidneys, liver, and pancreas. Creatine is then transported in blood to other organs, such as muscles and the brain. About 1% to 2% of muscle creatine is converted to creatinine daily. Because the amount of endogenous creatinine produced is proportional to muscle mass, the production varies with age and sex. Given the same glomerular filtration rate, men have a higher creatinine concentration than women, muscular individuals have higher levels than less muscular individuals, and younger persons higher levels than the elderly. Depending on the individual's meat intake, diet may influence the value by about 10%. On the whole, however, dietary fluctuation of creatinine intake cause only minor variation in daily creatinine excretion of the same person. High creatinine values are found in acute and chronic kidney insufficiency and dehydration.

Both chemical and enzymatic methods are used to measure creatinine in body fluids. Under routine conditions, creatinine can be enzymatically determined with only a few interferences.

PRINCIPLE OF THE PROCEDURE

Creatinine is determined by enzymatic colorimetric assay as follows:

Creatinine is converted to sarcosine with the aid of creatininase and creatinase. Sarcosine is then converted to glycine, formaldehyde and hydrogen peroxide in the presence of oxygen by sarcosine oxidase. The liberated hydrogen peroxide reacts with 4-aminophenazone and HTIB to form a quinone imine chromogen in a reaction catalyzed by peroxidase. The color intensity is directly proportional to the concentration of creatinine present and can be measured photometrically at 540nm.

REAGENT INFORMATION

Reagent A 4 x 40 ml
Reagent B 4 x 20 ml

Concentrations

Reagent A	
TAPS** buffer, pH 8.1	30 mmol/l
Creatininase (microorganisms)	> 333 µkat/l
Sarcosine oxidase (microorganisms)	> 133 µkat/l
Ascorbate oxidase (microorganisms)	> 33 µkat/l
HTIB	5.9 mmol/l
Detergents	
Preservative	
Reagent B:	
TAPS** buffer, pH 8.0	50 mmol/l
Creatininase (microorganisms)	> 500 µkat/l
Peroxidase (horseradish)	> 16.7 µkat/l
4-aminophenazone	2.0 mmol/l
Potassium hexacyanoferrate (II)	163 µmol/l
Detergent	
Preservative	

** TAPS = 3-[[tris(hydroxymethyl)methyl]amino]propanesulfonic acid

Precautions

For *in vitro* diagnostic use only. Exercise the normal precautions required for handling all laboratory reagents.

Preparation

The reagents are ready for use.

Note: Check that there are no bubbles in the bottleneck or on the surface of the reagent when you insert the reagent vials or vessels in the Konelab analyzer.

Storage and Stability

Reagents in unopened vials are stable at 2...8 °C until the expiry date printed on the label. Refer to the Application Notes of your Konelab analyzer for the on board stability of reagents.

SPECIMEN COLLECTION**Sample Type**

Only the specimens given below were tested and found acceptable:
Serum, Li-heparin plasma, Na-EDTA plasma and urine.

NOTE for specimen collection and preparation, use only suitable tubes or collection containers. Collect urine without additives.

Precautions

Human samples should be handled and disposed of as if they were potentially infectious.

Storage (3)

The serum and plasma samples can be stored for 7 days at 20...25 °C or at 4...8 °C or for 3 months at -20 °C.

The urine sample can be stored for 2 days at 20...25 °C, for 6 days at 4...8 °C or for 6 months at -20 °C.

TEST PROCEDURE

Refer to the Reference Manual and Application Notes separate for serum/plasma and urine for an automated procedure on your Konelab analyzer. See the latest applications from the web page. Any application which has not been validated by Thermo Fisher Scientific Oy cannot be performance guaranteed and therefore must be evaluated by the user.

Materials provided

Reagents as described above.

Materials required but not provided

WashFluid, code: 981842.
Calibrator and controls as indicated below.

Calibration

Use sCal, code 981831 according to the instructions provided for your Konelab analyzer. Allow the calibrator to remain on the analyzer for a maximum of one hour.

Traceability:

Refer to the package insert of sCal.
The response (A) is converted to result units by a calculation factor.

Quality Control

Use quality control samples at least once a day and after each calibration and every time a new bottle of reagent is used. It is recommended to use two level controls.

Available controls:

Serum/plasma:
Nortrol, code: 981043
Abtrol, code: 981044


Urine:

uTrol, code: 981821
uTrol High, code: 981822

The Control intervals and limits must be adapted to the individual laboratory requirements. The results of the quality control samples should fall within the limits pre-set by the laboratory.

CALCULATION OF RESULTS

The results are calculated automatically by the Konelab analyzer using a calibration curve.

Calibration Curve (example)

Konelab 20/20XT/30/60. The calibration curve is lot dependent.

LIMITATIONS OF THE PROCEDURE**Interference**

Criterion: Recovery within ± 10% of initial values.

Serum/plasma:

Unconjugated Bilirubin: No interference found up to 400 µmol/l (23.4 mg/dl). Higher bilirubine concentrations cause erroneously low creatinine values.

Conjugated Bilirubin: No interference found up to 300 µmol/l (17.5 mg/dl). Higher bilirubine concentrations cause erroneously low creatinine values.

Hemoglobin / Hemolysis: No interference found up to 10 g/l of hemoglobin in hemolysate.

Lipemia: No interference found up to 10 g/l of Intralipid™ (trademark of Fresenius Kabi AB) or 12 mmol/l (1062mg/dl) of triglycerides. There is a poor correlation between turbidity and concentration of triglycerides.

Monoclonal immunoglobulins (or parts of them) may interfere with the assay. Results from patiens suspected of having such antibodies should be carefully evaluated.

For other interfering substances, please refer to the reference 4.

Urine:

Conjugated Bilirubin: No interference found up to 1000 µmol/l (58 mg/dl).

Hemoglobin / Hemolysis: No interference found up to 10 g/l of hemoglobin in hemolysate.

Glucose: No interference found up to 139 mmol/l (2500 mg/dl).

Ascorbic acid: No interference found up to 5.7 mmol/l (100 mg/dl).

Calcium dobesilate and α-methyl dopa cause erroneously low creatinine values.

For other interfering substances, please refer to the reference 4.

EXPECTED VALUES (5)**Serum/plasma:**

Male: 59 – 104 µmol/l (0.67 – 1.17 mg/dl)

Female: 45 – 84 µmol/l (0.51 – 0.95 mg/dl)

Urine (1st morning urine):

Male: 3.54 – 24.6 mmol/l (40 – 278 mg/dl)

Female: 2.55 – 20.0 mmol/l (29 – 226 mg/dl)

The quoted values should serve as a guide only. It is recommended that each laboratory verify this range or derive a reference interval for the population that it serves.

Conversion factor:

µmol/l x 0.0113 → mg/dl

mmol/l x 11.3 → mg/dl

MEASURING RANGE**Serum/plasma:**

10 – 2500 µmol/l (0.11 – 28 mg/dl)

Extended measuring range after secondary dilution: 10 – 10000 µmol/l (0.11 – 113 mg/dl)

Urine:

0.2 – 40 mmol/l (2.3 – 452 mg/dl)

Extended measuring range after secondary dilution: 0.01 – 200 mmol/l (0.11 – 2260 mg/dl)

PERFORMANCE CHARACTERISTICS

The results obtained in individual laboratories may differ from the performance data given.

Detection limit**Serum/plasma:**

2 µmol/l (0,02 mg/dl)

The detection limit represents the lowest measurable concentration/activity that can be distinguished from zero. It is calculated as the concentration of zero sample +3 SD (within run, n = 24).

Urine:

0,002 mmol/l (0,02 mg/dl)

The detection limit represents the lowest measurable concentration/activity that can be distinguished from zero. It is calculated as the concentration of zero sample +3 SD (within run, n = 24)

Imprecision**Serum/plasma:**

	Mean 38 µmol/l		Mean 155 µmol/l		Mean 484 µmol/l	
	SD	CV%	SD	CV%	SD	CV%
Within run	0,56	1,4	0,69	0,4	1,85	0,4
Between run	0,16	0,4	0,71	0,5	1,30	0,3
Total	0,84	2,1	2,32	1,5	6,76	1,4

A precision study was performed according to guidelines in the CLSI (former NCCLS) Document EP5-A varying calibrations and operators using one Konelab 60 for 20 days, with the number of measurements being n = 80.

Urine:

	Mean 6,7 mmol/l		Mean 14,6 mmol/l	
	SD	CV%	SD	CV%
Within run	0,06	1,0	0,13	0,9
Between run	0,05	0,8	0,17	1,2
Total	0,24	3,5	0,51	3,4

A precision study was performed according to guidelines in the CLSI (former NCCLS) Document EP5-A varying calibrations and operators using one Konelab 60 for 20 days, with the number of measurements being n = 80.

Method comparison**Serum/plasma:**

A comparison study was performed according to guidelines in the CLSI (former NCCLS) Document EP9-A using a commercially available enzymatic method as a reference.

Linear regression (result unit µmol/l):

$$y = 0,97 x - 0,3$$

$$r = 0,999$$

$$n = 110$$

The sample concentrations were between 38 and 1951 µmol/l.

Urine:

A comparison study was performed according to guidelines in the CLSI (former NCCLS) Document EP9-A using a commercially available enzymatic method as a reference.

Linear regression (result unit mmol/l):

$$y = 1,10 x + 0,02$$

$$r = 0,999$$

$$n = 119$$

The sample concentrations were between 1,4 and 33,2 mmol/l.

BIBLIOGRAPHY

- Burtis, CA and Ashwood, E R (ed.), Tietz Fundamentals of Clinical Chemistry, 5th edition, W B Saunders Company, Philadelphia, 2001, pp. 419-422.
- Thomas, L. (ed.), Clinical Laboratory Diagnostics; Use and Assessment of Clinical Laboratory Results, 1st edition, TH-Books Verlagsgesellschaft mbH, Frankfurt/Main, Germany, 1998, pp. 366-370.
- Guder WG, Narayanan S, Wisser H, Zawta B. List of Analytes; Preanalytical variables. Brochure in: Samples: From Patient to the Laboratory. GIT Verlag GmbH, Darmstadt, 1996.
- Young, D.S., Effects of Drugs on Clinical Laboratory Tests, Fifth Edition, AACCC Press, Washington, D.C., 2000, pp. 3-240 - 3-261.
- Mazzachi BC, Peake MJ, Ehrhardt V. Reference range and method comparison studies for enzymatic and Jaffe creatinine assays in plasma and serum and early morning urine. Clin Lab. 2000; 46: 53-55.

MANUFACTURED BY

Thermo Fisher Scientific Oy
Clinical Diagnostics Finland
Ratastie 2, P.O. Box 100, FI-01621 Vantaa, Finland
Tel. +358 9 329 100, Fax +358 9 3291 0300
www.thermo.com/konelab

Date of revision (yyyy-mm-dd)

2007-08-07

Changes from previous version

Company name updated.

DE

Konelab™ / T Series CREATININE (Enzymatic)

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BIBLIOGRAPHIE

- 1) Burtis, CA and Ashwood, E R (ed.), Tietz Fundamentals of Clinical Chemistry, 5th edition, W B Saunders Company, Philadelphia, 2001, pp. 419-422.
- 2) Thomas, L. (ed.), Clinical Laboratory Diagnostics; Use and Assessment of Clinical Laboratory Results, 1st edition, TH-Books Verlagsgesellschaft mbH, Frankfurt/Main, Germany, 1998, pp. 366-370.
- 3) Guder WG, Narayanan S, Wisser H, Zawta B. List of Analytes; Preanalytical variables. Brochure in : Samples: From Patient to the Laboratory. GIT Verlag GmbH, Darmstadt, 1996.
- 4) Young, DS. Effects of Drugs on Clinical Laboratory Tests, cinquième édition, AACCC Press, Washington, D.C., 2000, pp. 3-240 - 3-261.
- 5) Mazzachi BC, Peake MJ, Ehrhardt V. Reference range and method comparison studies for enzymatic and Jaffé creatinine assays in plasma and serum and early morning urine. Clin Lab. 2000; 46: 53-55.

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SEZNAM LITERATURY

- Burtis, CA and Ashwood, E R (ed.), Tietz Fundamentals of Clinical Chemistry, 5th edition, W Saunders Company, Philadelphia, 2001, pp. 419-422.
- Thomas, L. (ed.), Clinical Laboratory Diagnostics; Use and Assessment of Clinical Laboratory Results, 1st edition, TH-Books Verlagsgesellschaft mbH, Frankfurt/Main, Germany, 1998, pp. 366-370.
- Guder WG, Narayanan S, Wisser H, Zawta B. List of Analytes; Preanalytical variables. Brochure in: Samples: From Patient to the Laboratory. GIT Verlag GmbH, Darmstadt, 1996.
- Young, D.S., Effects of Drugs on Clinical Laboratory Tests, Fifth Edition, AACCC Press, Washington, D.C., 2000, pp. 3-240 - 3-261.
- Mazzachi BC, Peake MJ, Ehrhardt V. Reference range and method comparison studies for enzymatic and Jaffe creatinine assays in plasma and serum and early morning urine. Clin Lab. 2000; 46: 53-55.

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- Burtis, CA and Ashwood, E R (ed.), Tietz Fundamentals of Clinical Chemistry, 5th edition, W B Saunders Company, Philadelphia, 2001, pp. 419-422.
- Thomas, L. (ed.), Clinical Laboratory Diagnostics; Use and Assessment of Clinical Laboratory Results, 1st edition, TH-Books Verlagsgesellschaft mbH, Frankfurt/Main, Germany, 1998, pp. 366-370.
- Guder WG, Narayanan S, Wissner H, Zawta B. List of Analytes; Preanalytical variables. Brochure in: Samples: From Patient to the Laboratory. GIT Verlag GmbH, Darmstadt, 1996.
- Young, D.S., Effects of Drugs on Clinical Laboratory Tests, Fifth Edition, AACCC Press, Washington, D.C., 2000, pp. 3-240 - 3-261.
- Mazzachi BC, Peake MJ, Ehrhardt V. Reference range and method comparison studies for enzymatic and Jaffe creatinine assays in plasma and serum and early morning urine. Clin Lab. 2000; 46: 53-55.

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KIRJANDUS

- Burtis, CA and Ashwood, E R (ed.), Tietz Fundamentals of Clinical Chemistry, 5th edition, W B Saunders Company, Philadelphia, 2001, pp. 419-422.
- Thomas, L. (ed.), Clinical Laboratory Diagnostics; Use and Assessment of Clinical Laboratory Results, 1st edition, TH-Books Verlagsgesellschaft mbH, Frankfurt/Main, Germany, 1998, pp. 366-370.
- Guder WG, Narayanan S, Wissner H, Zawta B. List of Analytes; Preanalytical variables. Brochure in: Samples: From Patient to the Laboratory, GIT Verlag GmbH, Darmstadt, 1996.
- Young, D.S., Effects of Drugs on Clinical Laboratory Tests, Fifth Edition, AACC Press, Washington, D.C., 2000, pp. 3-240 - 3-261.
- Mazzachi BC, Peake MJ, Ehrhardt V. Reference range and method comparison studies for enzymatic and Jaffé creatinine assays in plasma and serum and early morning urine. Clin Lab. 2000; 46: 53-55.

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BIBLIOGRÁFIA

- Burtis, CA and Ashwood, E R (ed.), Tietz Fundamentals of Clinical Chemistry, 5th edition, W B Saunders Company, Philadelphia, 2001, pp. 419-422.
- Thomas, L. (ed.), Clinical Laboratory Diagnostics; Use and Assessment of Clinical Laboratory Results, 1st edition, TH-Books Verlagsgesellschaft mbH, Frankfurt/Main, Germany, 1998, pp. 366-370.
- Guder WG, Narayanan S, Wissner H, Zawta B. List of Analytes; Preanalytical variables. Brochure in: Samples: From Patient to the Laboratory, GIT Verlag GmbH, Darmstadt, 1996.
- Young, D.S., Effects of Drugs on Clinical Laboratory Tests, Fifth Edition, AACC Press, Washington, D.C., 2000, pp. 3-240 - 3-261.
- Mazzachi BC, Peake MJ, Ehrhardt V. Reference range and method comparison studies for enzymatic and Jaffé creatinine assays in plasma and serum and early morning urine. Clin Lab, 2000; 46: 53-55.

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BIBLIOGRAFIA

- Burtis, CA and Ashwood, E R (ed.), Tietz Fundamentals of Clinical Chemistry, 5th edition, W B Saunders Company, Philadelphia, 2001, pp. 419-422.
- Thomas, L. (ed.), Clinical Laboratory Diagnostics: Use and Assessment of Clinical Laboratory Results, 1st edition, TH-Books Verlagsgesellschaft mbH, Frankfurt/Main, Germany, 1998, pp. 366-370.
- Guder WG, Narayanan S, Wissler H, Zawta B. List of Analytes; Preanalytical variables. Brochure in: Samples: From Patient to the Laboratory. GIT Verlag GmbH, Darmstadt, 1996.
- Young, D.S., Effects of Drugs on Clinical Laboratory Tests, Fifth Edition, AACCPress, Washington, D.C., 2000, pp. 3-240 - 3-261.
- Mazzachi BC, Peake MJ, Ehrhardt V. Reference range and method comparison studies for enzymatic and Jaffe creatinine assays in plasma and serum and early morning urine. Clin Lab. 2000; 46: 53-55.

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- Burtis, CA and Ashwood, E R (ed.), Tietz Fundamentals of Clinical Chemistry, 5th edition, W B Saunders Company, Philadelphia, 2001, pp. 419-422.

- Thomas, L. (ed.), Clinical Laboratory Diagnostics; Use and Assessment of Clinical Laboratory Results, 1st edition, TH-Books Verlagsgesellschaft mbH, Frankfurt/Main, Germany, 1998, lpp. 366-370.
- Guder WG, Narayanan S, Wisser H, Zawta B. List of Analytes; Preanalytical variables. Brochure in: Samples: >From Patient to the Laboratory. GIT Verlag GmbH, Darmstadt, 1996.
- Young, D.S., Effects of Drugs on Clinical Laboratory Tests, Fifth Edition, AACCPress, Washington, D.C., 2000, pp. 3-240 - 3-261.
- Mazzachi BC, Peake MJ, Ehrhardt V. Reference range and method comparison studies for enzymatic and Jaffe creatinine assays in plasma and serum and early morning urine. Clin Lab. 2000; 46: 53-55.

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- 1) Burtis, CA and Ashwood, E R (ed.), Tietz Fundamentals of Clinical Chemistry, 5th edition, W B Saunders Company, Philadelphia, 2001, pp. 419-422.
- 2) Thomas, L. (ed.), Clinical Laboratory Diagnostics; Use and Assessment of Clinical Laboratory Results, 1st edition, TH-Books Verlagsgesellschaft mbH, Frankfurt/Main, Germany, 1998, pp. 366-370.
- 3) Guder WG, Narayanan S, Wissler H, Zawta B. List of Analytes; Preanalytical variables. Brochure in: Samples: From Patient to the Laboratory. GIT Verlag GmbH, Darmstadt, 1996.
- 4) Young, D.S., Effects of Drugs on Clinical Laboratory Tests, Fifth Edition, AACC Press, Washington, D.C., 2000, pp. 3-240 - 3-281.
- 5) Mazzachi BC, Peake MJ, Ehrhardt V. Reference range and method comparison studies for enzymatic and Jaffe creatinine assays in plasma and serum and early morning urine. Clin Lab. 2000; 46: 53-55.

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LITERATURA

- Burtis, CA and Ashwood, E R (ed.), Tietz Fundamentals of Clinical Chemistry, 5th edition, W B Saunders Company, Philadelphia, 2001, pp. 419-422.
- Thomas, L. (ed.), Clinical Laboratory Diagnostics; Use and Assessment of Clinical Laboratory Results, 1st edition, TH-Books Verlagsgesellschaft mbH, Frankfurt/Main, Germany, 1998, pp. 366-370.
- Guder WG, Narayanan S, Wisser H, Zawta B. List of Analytes; Preanalytical variables. Brochure in: Samples: From Patient to the Laboratory. GIT Verlag GmbH, Darmstadt, 1996.
- Young, D.S., Effects of Drugs on Clinical Laboratory Tests, Fifth Edition, AACC Press, Washington, D.C., 2000, pp. 3-240 - 3-261.
- Mazzachi BC, Peake MJ, Ehrhardt V. Reference range and method comparison studies for enzymatic and Jaffe creatinine assays in plasma and serum and early morning urine. Clin Lab. 2000; 46: 53-55.

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BIBLIOGRAFIA

- Burtis, CA and Ashwood, E R (ed.), Tietz Fundamentals of Clinical Chemistry, 5th edition, W B Saunders Company, Philadelphia, 2001, pp. 419-422.
- Thomas, L. (ed.), Clinical Laboratory Diagnostics; Use and Assessment of Clinical Laboratory Results, 1st edition, TH-Books Verlagsgesellschaft mbH, Frankfurt/Main, Germany, 1998, pp. 366-370.
- Guder WG, Narayanan S, Wisser H, Zawta B. List of Analytes; Preanalytical variables. Brochure in: Samples: From Patient to the Laboratory. GIT Verlag GmbH, Darmstadt, 1996.
- Young, D.S., Effects of Drugs on Clinical Laboratory Tests, Fifth Edition, AAC Press, Washington, D.C., 2000, pp. 3-240 - 3-261.
- Mazzachi BC, Peake MJ, Ehrhardt V. Reference range and method comparison studies for enzymatic and Jaffé creatinine assays in plasma and serum and early morning urine. Clin Lab. 2000; 46: 53-55.

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REFERENSER

- 1) Burtis, CA and Ashwood, E R (ed.), Tietz Fundamentals of Clinical Chemistry, 5th edition, W B Saunders Company, Philadelphia, 2001, pp. 419-422.
- 2) Thomas, L. (ed.), Clinical Laboratory Diagnostics; Use and Assessment of Clinical Laboratory Results, 1st edition, TH-Books Verlagsgesellschaft mbH, Frankfurt/Main, Germany, 1998, pp. 366-370.
- 3) Guder WG, Narayanan S, Wisser H, Zawta B. List of Analytes; Preanalytical variables. Brochure in: Samples: From Patient to the Laboratory. GIT Verlag GmbH, Darmstadt, 1996.
- 4) Young, D.S., Effects of Drugs on Clinical Laboratory Tests, Fifth Edition, AACC Press, Washington, D.C., 2000, pp. 3-240 - 3-261.
- 5) Mazzachi BC, Peake MJ, Ehrhardt V. Reference range and method comparison studies for enzymatic and Jaffé creatinine assays in plasma and serum and early morning urine. Clin Lab. 2000; 46: 53-55.

TILLVERKARE

Thermo Fisher Scientific Oy
Clinical Diagnostics Finland
Ratastie 2, P.O. Box 100, FI-01621 Vantaa, Finland
Tel. +358 9 329 100, Fax +358 9 3291 0300
www.thermo.com/konelab

Revisionsdatum (åååå-mm-dd)

2007-08-07

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