



Purchase 1 case of NUNC™ Chamber Slides; get 1 FREE case of NUNC EZ Flip™ Tubes

Purchase one case of NUNC Lab-Tek™ II or Lab-Tek II CC² Chamber Slide™ System and receive a free case of 500 NUNC EZ Flip 15mL Conical Centrifuge Tubes (Cat. No. 12-565-802)



Why use Lab-Tek Chamber Slides?

Inoculate, incubate, fix and stain in one complete system

What is the difference between Lab-Tek II and Lab-Tek II CC²?

Lab-Tek II has a Soda Lime Glass slide that has been washed to promote cell growth and adherence.

Lab-Tek II CC² has a chemically modified glass growth surface with a positive charge that mimics Poly-D-Lysine.



- Removable, polystyrene medium chamber of 1-, 2-, 4-, or 8-well configuration
- Nonfluorescent microscope slide, glass (25 x 75 x 1.2mm) with rounded corners
- Biocompatible adhesive
- Polystyrene cover
- Inert hydrophobic well border printed on slide
- Superfrost™ printed writing area
- Slide separator in each pack
- CE marked

NUNC EZ Flip 15mL Conical Centrifuge Tubes

- Patented hinged-cap system
- Ergonomic design for one-handed opening and closing
- Integrated leakproof cap prevents cross contamination and lost closures
- Clarified polypropylene combines chemical resistance and high speed tolerance
- White writing area and flat caps for sample identification
- Printed graduations
- Standard dimensions compatible with centrifuge rotors

Cat. No.	Description	No. of Wells	Culture Area, cm ² /well	Case Quantity	Case List
12-565-5	LabTek II Chamber Slide	1	8.6	96	\$496.50
12-565-6	LabTek II Chamber Slide	2	4.0	96	\$524.00
12-565-7	LabTek II Chamber Slide	4	1.7	96	\$538.25
12-565-8	LabTek II Chamber Slide	8	0.7	96	\$554.50
12-565-4	LabTek II CC ₂ Chamber Slide	1	8.6	96	\$730.50
12-565-3	LabTek II CC ₂ Chamber Slide	2	4.0	96	\$755.00
12-565-2	LabTek II CC ₂ Chamber Slide	4	1.7	96	\$779.25
12-565-1	LabTek II CC ₂ Chamber Slide	8	0.7	96	\$804.75

To claim your free product, complete the below information and fax it with your invoice to Thermo Fisher Scientific at **585-586-3294**.

Name: _____

Institution/Company: _____

Business Address (No PO Boxes): _____

Bldg./Room/Lab: _____

City: _____ State: _____ Zip Code: _____

Telephone: _____

E-Mail Address: _____

Offer valid on purchases between July 1 and December 31, 2007. Submissions must be received by January 15, 2008. Limit five free cell scrapers per lab. Void where prohibited by law or company policy. Free product specifications are subject to change. Not valid on prior purchases.



For customer service, call 1-800-766-7000.
To fax an order, use 1-800-926-1166.
To order online: www.fishersci.com