


## Fisherbrand

Polypropylene Micro Test Tube Racks  
Interlocking Four-Way Tube Rack  
Microcentrifuge Tube Racks  
96-Well Rack for PCR Tubes  
Stainless-Steel Chest Freezer Racks  
Premium Microcentrifuge Tubes  
96-Well Microskirted PCR Plates  
96-Well PCR Tubes  
96-Well PCR Tubes

SureOne\* Pipet Tips  
Polypropylene Biohazard Autoclave Bags  
Cell Strainers  
Cell Lifter  
Cell Spreaders  
Cell Scrapers  
Polyurethane Ice Buckets  
Holders for Biohazard Autoclave Bags  
Electroporation Cuvettes Plus\*  
Fine-Tip Marking Pens  
Fisher Scientific Traceable\* Countdown Timers  
Inoculating Loops  
Raised Ridge Petri Dishes with Clear Lids  
Disposable Borosilicate Glass Tubes with Plain End  
Disposable Pestle System  
Disposable 50mL Pipet Basins  
Absorbent Surface Protector



For a comprehensive listing of genomics products, visit  
[www.fishersci.com/genomics](http://www.fishersci.com/genomics).

**Fisherbrand\* Polypropylene Micro Test Tube Racks**

**Autoclavable racks hold up to 20 tubes**

- Racks permit storage of 1.5mL micro test tubes in prenumbered cavities
- Special step design holds 10 tubes in top row and 10 tubes in bottom row



| Cat. No.  | Each    | Case of 6 |
|-----------|---------|-----------|
| 14-754-30 | \$18.44 | \$93.72   |

**Fisherbrand\* Interlocking Four-Way Tube Racks**

**Each autoclavable rack holds four 50mL conical tubes, twelve 15mL conical tubes, or thirty-two 0.5 and 1.5mL microtubes**

- Rugged polypropylene racks accommodate tubes from 0.5 to 50mL by simply turning the rack
- Easily link together for custom configurations



| Color | Cat. No.  | Pack of 5 |
|-------|-----------|-----------|
| Blue  | 03-448-12 | \$40.15   |

Also comes in natural, assorted, and fluorescent green, orange, pink and yellow

**Fisherbrand\* Premium Microcentrifuge Tubes**



**Flat-top snap caps provide a safe, liquid-tight seal even with prolonged boiling**

- High-clarity polypropylene tubes withstand forces up to 30,000xG, fit all standard rotors
- Easy-open caps have needle insertion spot and textured marking area


| Color   | Capacity | O.D. x L      | Cat. No.   | Pack of 500 |
|---------|----------|---------------|------------|-------------|
| Natural | 0.5mL    | 7.5 x 30.5mm  | 05-408-120 | \$51.12     |
|         | 1.5mL    | 10.8 x 40.6mm | 05-408-129 | \$52.60     |
|         | 2mL      | 11 x 40.6mm   | 05-408-138 | \$53.23     |

Also comes in blue, green, red, yellow, amber, orange, lavender and assortment

**Fisherbrand\* Microcentrifuge Tube Racks**

**Two frosted writing areas and alphanumeric index for sample identification**

- Reversible, autoclavable polypropylene rack withstands freezer storage
- Holds ninety-six 1.5mL microcentrifuge tubes on one side and screw-cap microcentrifuge tubes on the other



| Color | Cat. No.  | Pack of 5 | Case of 4 |
|-------|-----------|-----------|-----------|
| Black | 05-541-37 | \$75.04   | \$265.63  |

Also comes in natural, lavender, blue, green, orange, pink, yellow, nonfluorescent as-sorted and fluorescent green, orange and yellow

**Fisherbrand\* 96-Well Rack for PCR Tubes**

**Autoclavable polypropylene rack for 0.2mL tubes**

- Imprinted alphanumeric grid and removable lid
- Each pack contains five racks in assorted colors (blue, green, pink, yellow, and orange)
- Accommodates individual tubes or 8- or 12-tube strips in a 96-well array




| Cat. No.  | Pack of 5 | Case of 40 Pk. |
|-----------|-----------|----------------|
| 03-448-20 | \$41.25   | \$1,319.80     |

**Fisherbrand\* 96-Well Nonskirted PCR Plates**

**Compatible with all 0.2mL thermal cyclers and ABI PRISM\* 310, 3100, and 3700 capillary sequencers**

- Virgin polypropylene with clear well bottom, raised rim
- Cut-away corner and alphanumeric grid
- 300µL maximum volume when used with adhesive and heat seals



| Color   | Cat. No.   | Case of 25 |
|---------|------------|------------|
| Natural | 14-230-232 | \$84.89    |
| Blue    | 14-230-234 | \$88.40    |
| Green   | 14-230-235 | \$88.40    |
| Red     | 14-230-233 | \$88.40    |
| Yellow  | 14-230-236 | \$88.40    |

**Fisherbrand® Stainless-Steel Chest Freezer Racks**



**Organize and optimize the freezers in your lab**

- Store a variety of sample containers, including boxes, plates, and tubes
- Sturdy, rust-resistant stainless steel

| No. of Shelves/<br>Capacity<br>for 2" (5cm)<br>Boxes | L x H x D                               | Cat. No.    | Each     |
|--|---|-------------|----------|
| 4  | 8.9 x 5.6 x 5.5" (22.6 x 14.2 x 14cm)   | 03-395-632  | \$30.83  |
| 6  | 13.9 x 5.6 x 5.5" (78 x 14.2 x 14cm)    | 03-395-634  | \$45.20  |
| 9  | 19.8 x 5.6 x 5.5" (50.3 x 14.2 x 14cm)  | 03-395-631  | \$65.68  |
| 10   | 22.7 x 5.6 x 5.5" (57.6 x 14.2 x 14cm)  | 03-395-613† | \$103.62 |
| 11   | 24.25 x 5.6 x 5.5" (61.6 x 14.2 x 14cm) | 03-395-614  | \$109.40 |
| 12   | 26.4 x 5.6 x 5.5" (67 x 14.2 x 14cm)    | 03-395-615  | \$115.41 |
| 13   | 28.6 x 5.6 x 5.5" (72.6 x 14.2 x 14cm)  | 03-395-616  | \$121.95 |

†With handle.

**Fisherbrand® 0.2mL PCR Tube Strips**



**Ideal for use in 0.2mL, 96-well V-bottom thermal cyclers**

- Virgin polypropylene with thin-wall design
- Strips of eight tubes can be cut into individual tubes
- 0.25mL maximum tube volume when capped

| Color   | Style              | Cat. No.   | Case of 250 |
|---------|--------------------|------------|-------------|
| Natural | Flat Cap Strips    | 14-230-215 | \$401.39    |
|         | Domed Cap Strips   | 14-230-210 | \$402.54    |
|         | Without Cap Strips | 14-230-220 | \$380.44    |

Also comes in Red, Blue, Green and Yellow

**Fisherbrand® 0.5mL Flat-Cap PCR Tubes**

**With integral quotation "snap shut" cap**

- Virgin polypropylene with thin-wall design
- Suitable for 0.5mL thermal cycler blocks
- 0.6mL maximum volume when closed



| Color   | Cat. No.   | Case of 1000 |
|---------|------------|--------------|
| Natural | 14-230-200 | \$148.25     |

Also comes in Red, Blue, Green and Yellow

**Fisherbrand® 0.2mL Flat- and Domed-Cap PCR Tubes**

**With integral quotation "snap shut" cap**


- Virgin polypropylene with thin-wall design
- Suitable for 0.2mL thermal cycler blocks
- 0.25mL maximum volume when closed



| Color   | Style           | Cat. No.   | Case of 1000 |
|---------|-----------------|------------|--------------|
| Natural | With Flat Caps  | 14-230-225 | \$154.84     |
|         | With Domed Caps | 14-230-205 | \$111.89     |

Also comes in Red, Blue, Green and Yellow

**Fisherbrand® Flat and Domed 8-Cap Strips**



**Suitable for use with 96-well PCR plates and PCR tube strips**

- Virgin polypropylene
- Eight caps per strip
- Flat caps offer large, labeled end tabs for easy removal and strip orientation

| Color   | Style     | Cat. No.   | Case of 250 |
|---------|-----------|------------|-------------|
| Natural | Flat Cap  | 14-230-230 | \$118.35    |
|         | Domed Cap | 14-230-231 | \$119.53    |

**Fisherbrand® Adhesive Plate Seals**

**Seals are optically clear for PCR applications**

| Cat. No.   | Case of 100 |
|------------|-------------|
| 08-408-240 | \$94.68     |



**Fisherbrand\* 96-Well Elastomer Cap Mat**

For 1 and 2mL DeepWell\* polypropylene microplates

- Color: Natural
- Pierceable, nonsterile
- Can be used at temperatures down to -80°C

| Cat. No.   | Case of 50 |
|------------|------------|
| 12-566-122 | \$228.89   |

**Fisherbrand\* MicroMat Cap Mat**

For pressure-sealing thermal cyclers

- For use with 96-well plates

| Cat. No.   | Pack of 5 |
|------------|-----------|
| 08-408-219 | \$51.00   |



**Fisherbrand\* Silicone Plate Seals**

Excellent chemical compatibility and resistance to tearing and coring

- Nonslit seals reseal after multiple penetrations to prevent sample evaporation and cross-well contamination
- Apply by hand or with an applicator
- Remain flexible at ambient temperatures and in the absence of chlorinated solvents



| Fits   | Cat. No.  | Pack of 10 |
|--|-----------|------------|
| Costar* 96-Well Polypropylene Plates and Assay Blocks, Whatman* 96-Well Multi-Chem* Microplates                  | 03-396-49 | \$101.40   |
| Fisherbrand Round-Bottom 96-Well Polypropylene Plates and Assay Blocks, Nunc* 96-Well Storage Plates             | 03-396-51 | \$101.40   |
| Fisherbrand 96-Well 2.2mL Assay Blocks, Whatman 96-Well Polypropylene Storage Plates (Square-Well Plate)         | 03-396-53 | \$109.72   |
| Fisherbrand 384-Well Polypropylene Plates, Costar Brand 384-Well Polypropylene Storage Block (Square-Well Plate) | 03-396-55 | \$116.45   |

**Fisherbrand\* 81-place Cryo Boxes**

Taller than other available storage boxes to accommodate vials or tubes up to 2.13" (5.35cm) high

- Hold 0.5, 1.5, and 2.0mL microtubes as well as most brands of vials
- Stackable; clear bottom; alphanumeric grid inside and outside
- Withstand conditions from -112°F to +249.8°F (-80°C to +121°C)



| Lid Color | Cat. No.   | Pack of 5 | Case of 10 Pk. |
|-----------|------------|-----------|----------------|
| Blue      | 03-395-402 | \$45.54   | \$389.22       |
| Green     | 14-826-340 | \$45.54   | \$389.22       |
| Natural   | 03-395-400 | \$45.54   | \$389.22       |
| Orange    | 03-395-404 | \$45.54   | \$389.22       |
| Purple    | 03-395-403 | \$45.54   | \$389.22       |

**Fisherbrand\* Polypropylene Biohazard Autoclave Bags**



Supplied with autoclavable elastic steam-venting closures

- Strong, pliable, leak- and puncture-resistant
- Clearly visible biohazard symbol and sterilization indicator
- Maximum temperature, 140°C (285°F); for use at 134°C (274°F)
- 2mil thick

| W x L                                   | Cat. No. | Pack of 200 |
|---|----------|-------------|
| <b>Printed Red Polypropylene Bags</b>   |          |             |
| 8 x 12" (20 x 30cm)                     | 01-828A  | \$109.87    |
| 14 x 19" (36 x 48cm)                    | 01-828B  | \$179.50    |
| 19 x 23" (48 x 58cm)                    | 01-828C  | \$243.36    |
| 25 x 35" (64 x 89cm)                    | 01-828D  | \$382.06    |
| 37 x 48" (94 x 122cm)                   | 01-828E  | \$371.06    |
| <b>Plain Red Polypropylene Bags</b>     |          |             |
| 8 x 12" (20 x 30cm)                     | 01-830A  | \$90.94     |
| 14 x 19" (36 x 48cm)                    | 01-830B  | \$165.47    |
| 19 x 23" (48 x 58cm)                    | 01-830C  | \$218.53    |
| 25 x 35" (64 x 89cm)                    | 01-830D  | \$332.50    |
| 37 x 48" (94 x 122cm)                   | 01-830E  | \$284.50    |
| <b>Printed Clear Polypropylene Bags</b> |          |             |
| 12 x 24" (30 x 61cm)                    | 01-826A  | \$187.39    |
| 24 x 30" (61 x 76cm)                    | 01-826B  | \$348.10    |
| 24 x 30" (61 x 76cm)                    | 01-826C  | \$405.74    |
| 31 x 38" (79 x 97cm)                    | 01-826D  | \$540.18    |
| 37 x 48" (94 x 122cm)                   | 01-826E  | \$394.37    |
| <b>Plain Clear Polypropylene Bags</b>   |          |             |
| 8 x 12" (20 x 30cm)                     | 01-826-1 | \$86.48     |
| 10 x 10" (25 x 25cm)                    | 01-826-2 | \$79.89     |
| 10 x 12" (25 x 30cm)                    | 01-826-3 | \$85.03     |
| 10 x 15" (25 x 38cm)                    | 01-826-4 | \$100.56    |
| 12 x 24" (30 x 61cm)                    | 01-826-5 | \$146.52    |
| 24 x 30" (61 x 76cm)                    | 01-826-6 | \$330.74    |
| 24 x 30" (61 x 76cm)                    | 01-826-7 | \$353.51    |

**Fisherbrand\* Cryo/Freezer Boxes**



| L x W x H                                  | No. of Cells | Cat. No.   | Each   |
|--|--------------|------------|--------|
| <b>With Dividers</b>                       |              |            |        |
| 3 x 3 x 2" (7.6 x 7.6 x 5.1cm)             | 25           | 03-395-112 | \$6.50 |
| 5.25 x 5.25 x 2" (13.4 x 13.4 x 5.1cm)     | 81           | 03-395-113 | \$8.44 |
| 5.25 x 5.25 x 2" (13.4 x 13.4 x 5.1cm)     | 100          | 03-395-114 | \$8.44 |
| <b>Without Dividers</b>                    |              |            |        |
| 3 x 3 x 2" (7.6 x 7.6 x 5.1cm)             | ----         | 03-395-100 | \$3.25 |
| 5.25 x 5.25 x 2" (13.4 x 13.4 x 5.1cm)     | ----         | 03-395-102 | \$4.36 |
| 5.25 x 5.25 x 3" (13.4 x 13.4 x 7.6cm)     | ----         | 03-395-103 | \$6.51 |
| 5.75 x 5.75 x 4.88" (14.6 x 14.6 x 12.3cm) | ----         | 03-395-101 | \$7.57 |

**Special moisture-repellent coating protects samples from the effects of liquid nitrogen or condensation from refrigerated equipment**

- Can be used in mechanical freezers and liquid nitrogen tanks at temperatures to -196°C
- Protective coating prevents water or liquid nitrogen from saturating box for increased longevity
- Bottom drain holes allow safe drainage of liquid nitrogen, prevent moisture buildup

**Fisherbrand\* SureOne\* Pipet Tips**

**Both racked sterile and racked or bulk nonsterile nonfiltered pipet tips are certified free of RNase, DNase and DNA**

- Racked sterile tips are shrink wrapped and certified free of ATP, pyrogen, bioburden, PCR inhibitors and endotoxin
- Universal fit<sup>1</sup>



| Volume  | Length | Packaging                     | Color   | Cat. No.   | Pack of 960 | Case of     |
|---|--------|-------------------------------|---------|------------|-------------|-------------|
| <b>Micropoint Tip, Graduated at 2µL</b>                   |        |                               |         |            |             |             |
| 0.01 to 10µL  | 1.25"  | Hinged Rack, 10 x 96          | Natural | 02-707-436 | \$45.05     | 10/\$222.15 |
| 0.01 to 10µL  | 1.25"  | Hinged Rack, Sterile, 10 x 96 | Natural | 02-707-437 | \$49.50     | 5/\$245.35  |
| 0.01 to 10µL  | 1.25"  | Reload Decks, 10 x 96         | Natural | 02-707-438 | \$40.00     | 5/\$197.45  |
| <b>Micropoint Tip, Eppendorf Style</b>                    |        |                               |         |            |             |             |
| 0.5 to 10µL   | 1.25"  | Hinged Rack, 10 x 96          | Natural | 02-707-433 | \$45.05     | 10/\$222.15 |
| 0.5 to 10µL   | 1.25"  | Hinged Rack, Sterile, 10 x 96 | Natural | 02-707-434 | \$49.50     | 5/\$245.30  |
| <b>Beveled Tip</b>  |        |                               |         |            |             |             |
| 1 to 200µL  | 2"     | Hinged Rack, 10 x 96          | Natural | 02-707-418 | \$42.10     | 5/\$208.05  |
| 1 to 200µL  | 2"     | Hinged Rack, Sterile, 10 x 96 | Natural | 02-707-419 | \$47.10     | 5/\$233.25  |
| 1 to 200µL  | 2"     | Reload Decks, 10 x 96         | Natural | 02-707-420 | \$35.55     | 5/\$176.40  |
| 1 to 200µL  | 2"     | Hinged Rack, 10 x 96          | Yellow  | 02-707-421 | \$42.10     | 5/\$208.05  |
| 1 to 200µL  | 2"     | Hinged Rack, Sterile, 10 x 96 | Yellow  | 02-707-422 | \$47.10     | 5/\$233.25  |
| 1 to 200µL  | 2"     | Reload Decks, 10 x 96         | Yellow  | 02-707-423 | \$35.55     | 5/\$176.40  |
| <b>Micropoint Tip, Graduated at 100, 200, 500, 1000µL</b> |        |                               |         |            |             |             |
| 100 to 1250µL   | 2.75"  | Hinged Rack, 10 x 96          | Natural | 02-707-401 | \$54.00     | 5/\$267.20  |
| 100 to 1250µL   | 2.75"  | Hinged Rack, Sterile, 10 x 96 | Natural | 02-707-402 | \$62.15     | 5/\$306.40  |
| 100 to 1250µL   | 2.75"  | Rack Inserts, 10 x 96         | Natural | 02-707-403 | \$49.25     | 5/\$243.25  |
| 100 to 1250µL   | 2.75"  | Hinged Rack, 10 x 96          | Blue    | 02-707-406 | \$54.00     | 5/\$267.20  |
| 100 to 1250µL   | 2.75"  | Hinged Rack, Sterile, 10 x 96 | Blue    | 02-707-407 | \$62.15     | 5/\$306.40  |
| 100 to 1250µL   | 2.75"  | Rack Inserts, 10 x 96         | Blue    | 02-707-408 | \$49.25     | 5/\$243.25  |

<sup>1</sup>If more detailed information regarding a specific pipetter is required, please contact your Fisher Scientific Customer Service Representative. Additional sizes and packaging available

**Fisherbrand\* Cell Strainers**



**Universal fit for compatibility with all 50mL conical tubes**

- Provide a fast, simple alternative to gauze filtration when dissociating cells from clumps or primary tissues
- Produce more uniform single-cell suspensions
- Polypropylene frame with nylon mesh

| Mesh Size | Color  | Cat. No.   | Case of 50 |
|-----------|--------|------------|------------|
| 40µm      | Blue   | 22-363-547 | \$99.00    |
| 70µm      | White  | 22-363-548 | \$99.00    |
| 100µm     | Yellow | 22-363-549 | \$99.00    |

**Fisherbrand\* Cell Lifter†**

**Disposable lifter quickly removes cell layers undamaged**



- Maintain even contact with growth area in culture dishes
- Made entirely of polyethylene
- Presterilized and individually wrapped

| Handle | Blade | Cat. No. | Case of 100 |
|--------|-------|----------|-------------|
| 18cm   | 2cm   | 08-773-1 | \$162.51    |

**Fisherbrand\* Cell Scrapers**



**Gently harvest cells from tissue culture flasks**

- Polystyrene handles are compatible with a variety of culture flasks
- Polyethylene blades swivel for even contact with growth surface
- Presterilized and individually wrapped

| Handle Length | Blade Length | Cat. No. | Case of 100 |
|---------------|--------------|----------|-------------|
| 25cm          | 1.8cm        | 08-773-2 | \$235.54    |
| 39cm          | 3.0cm        | 08-773-3 | \$256.24    |

**Fisherbrand\* Cell Spreaders**



**Nonbreakable and reusable polypropylene spreaders available in two sizes: 30 and 60mm bar widths**

- Use 30mm wide spreaders for prepared biplates; 60mm for 100 to 150mm prepared plates
- Autoclavable

| Size              | Color  | Cat. No.  | Pack of 25 | Case of 4 Pack |
|-------------------|--------|-----------|------------|----------------|
| <b>Nonsterile</b> |        |           |            |                |
| 30mm              | Blue   | 08-100-10 | \$32.89    | \$102.89       |
| 60mm              | Yellow | 08-100-11 | \$29.64    | \$92.55        |

**Fisherbrand\* Polyurethane Ice Buckets**

**For use with liquid nitrogen, dry ice, water ice, and salt solutions**



- Rigid polyurethane skinned-foam construction provides high durability and excellent insulating properties; buckets will not sweat or leak, and are unbreakable with normal use
- Resistant to acetone and easy to clean—can be chemically disinfected
- Withstand temperatures up to 95°C (203°F)

| Capacity | Color | Cat. No.  | Each     |
|----------|-------|-----------|----------|
| 4.5L     | Black | 02-591-44 | \$103.00 |
| 4.5L     | Blue  | 02-591-45 | \$103.00 |
| 4.5L     | Red   | 02-591-46 | \$103.00 |



## Fisherbrand® Fine-Tip Marking Pens



Designed to write on cold and damp surfaces as easily as they write on dry surfaces. Proven permanent on most surfaces, including plastic, metal, and ceramic

- Quick-drying ink has no unpleasant odor
- Will not smear when subjected to water or alcohol during normal use
- Highly resistant to abrasion

| Ink Color             | Cat. No. | Pack of 10 |
|-----------------------|----------|------------|
| <b>Extra-Fine Tip</b> |          |            |
| Black                 | 13-379-1 | \$41.80    |

Also available in red, blue, green

## Fisherbrand® Holders for Biohazard Autoclave Bags



Hold Fisherbrand and most biohazard autoclave bags

- High-density polyethylene container with rigid low-density polyethylene lid with knob
- Lid fits securely to seal against hazardous aerosols
- Containers have integral handle for easy transporting

| Outside D x H           | Accepts Bag Size | Capacity | Cat. No.   | Each    | Case of 6 |
|-------------------------|------------------|----------|------------|---------|-----------|
| 9 x 11"<br>(23 x 28cm)  | 12" (30cm)       | 8 qt.    | 14-375-175 | \$31.75 | \$162.23  |
| 12 x 16"<br>(30 x 41cm) | 16" (40.6cm)     | 24 qt.   | 14-375-176 | \$54.75 | \$280.56  |
| 15 x 20"<br>(38 x 51cm) | 20" (50.8cm)     | 48 qt.   | 14-375-177 | \$67.59 | \$346.63  |

## Fisherbrand® Electroporation Cuvettes Plus\*

Disposable, color-coded cuvetts fit most electroporators. Round cap allows single-handed removal

- Clear, medical grade polycarbonate for high visibility
- Polished aluminum lowers arcing frequency
- Sterilized by gamma irradiation and supplied individually wrapped



| Description    | Cap Color | Cat. No. | Pack of 50 |
|----------------|-----------|----------|------------|
| 1mm gap, 90µL  | Purple    | FB101    | \$138.32   |
| 2mm gap, 400µL | Blue      | FB102    | \$138.32   |
| 4mm gap, 800µL | Yellow    | FB103    | \$138.32   |

## Fisher Scientific Traceable® Countdown Timers

For general-purpose timing applications

- One timing channel and a single-line display
- Rugged ABS plastic case
- Traceable to NIST\* standards




| Timer Range | LCD Display Height | Power Supply        | Alarm                | Cat. No.  | Each    |
|-------------|--------------------|---------------------|----------------------|-----------|---------|
| 100 hours   | 0.75"              | 2-year silver oxide | 1 min., high-decibel | 14-649-17 | \$30.68 |

## Fisherbrand® Inoculating Loops

No need to flame metal loops and expose to airborne pathogens

- Ideal for quantitative procedures such as sampling, serial dilutions and urine counts, as well as for bacterial inoculation
- Molded from biologically inert polymer
- Flexible enough to facilitate streaking, yet rigid enough to penetrate the media for subsurface inoculation



| Description           | Color      | Packaging | Cat. No.   | Pack of 960 | Case of 10 |
|-----------------------|------------|-----------|------------|-------------|------------|
| 10µL Loop; 3.8mm I.D. | Light Blue | 10/Pk.    | 22-363-597 | \$210.05    | \$1,912.45 |

Needles and Combi-loops also available.

**Fisherbrand\* Raised Ridge Petri Dishes with Clear Lids**



Polystyrene, radiation-sterilized dishes are easy to handle with just one hand. Stack high and remain stacked, reducing slipping

- Top and bottom edges are squared off (not rounded) to provide a firm, sure grip
- Bottom is flat so medium spreads evenly over the entire dish bottom

- Lid features three venting ribs, flared skirt; fits loosely for air exchange

| O.D. x H   | Cat. No.  | Case of 500 |
|------------|-----------|-------------|
| 100 x 15mm | 08-757-12 | \$179.00    |

**Fisherbrand\* Disposable Borosilicate Glass Tubes with Plain End**



Made of borosilicate glass to reduce pH changes and contaminants potentially leached from soda lime glass


- Ideal for tissue culture and clinical chemistry
- Packaged in shrink-wrapped trays of 250 or 125 to prevent contamination or breakage
- 10 and 12mm O.D. sizes are suitable for cell washing procedures

| O.D. x L  | Cat. No.  | Case of 1000 |
|-----------|-----------|--------------|
| 12 x 75mm | 14-961-26 | \$76.61      |

**Fisherbrand\* Disposable Pestle System**

Low-cost alternative to mechanical or metal and glass tissue grinders

- Lightweight, cordless mixer can be used with disposable pestles and microtubes
- High-speed vortexing action completely homogenizes cells and tissue within seconds; re-suspends protein and DNA pellets
- Rugged polybutylene terephthalate (PBTP) plastic pestle yields superior homogenization



| Cat. No.   | Each     |
|------------|----------|
| 03-392-106 | \$101.43 |

**Fisherbrand\* Disposable 50mL Pipet Basins**



Bulk packaged in easy-to-open "Snap Smartly" bags to reduce contamination risks

- Printed graduations for more accurate filling
- Easy-pour corners for quick disposal of excess sample solution
- Extended writing surface for product I.D.

| Material | Sterile | Cat. No.   | Case of     |
|----------|---------|------------|-------------|
| PVC      | No      | 13-681-500 | 100/\$71.30 |

**Fisherbrand\* Absorbent Surface Protector**

Disposable sheets of absorbent paper with polyethylene backing

- Ideal for preventing glass breakage, protecting benches, absorbing liquids, and lining humidity chambers and fume cupboards
- Paper side quickly absorbs liquids and keeps spills from spreading
- Polyethylene side prevents liquids from reaching work area and is useful for sample recovery



| Width x Length         | Cat. No.   | Pack of 50 | Case of 8 Pack |
|------------------------|------------|------------|----------------|
| 18 x 22.5" (46 x 57cm) | 12-007-187 | \$104.66   | \$754.69       |