

FALL FRENZY SPECIAL OFFER

Buy Any 2, Receive a 3rd* Free
Now Through Dec 30, 2005

Immobilon™ Western HRP Detection Reagent

- High sensitivity over a broad detection range
- Reduces antibody consumption by 2 to 5 fold
- Signal duration > 2 hours
- Working substrate stable for 10 days



Amicon® Ultra Centrifugal Devices

- High recovery Ultracel membrane in a range of molecular weight cutoffs
- Vertical membrane reduces concentration polarization for ultra-fast spin times of ~10 -15 minutes
- Patented deadstop for reproducible concentrate volumes



*Same or lesser value

Ordering Information

Immobilon Western HRP Detection Reagent *
 *AP substrate also available, please inquire

Amicon Ultra Centrifugal Devices 24/Pack*
 *96/pack also available, please inquire

		<u>4 mL Volume</u>		<u>15 mL volume</u>	
WBKLS0100	HRP Reagent (100 mL)	UFC800524	5 Kd	UFC900524	5 Kd
WBKLS0500	HRP Reagent (500 mL)	UFC801024	10 Kd	UFC901024	10 Kd
		UFC803024	30 Kd	UFC903024	30 Kd
		UFC805024	50 Kd	UFC905024	50 Kd
		UFC810024	100 Kd	UFC910024	100 Kd

To receive your free Detection Reagent or Amicon Ultra, simply fill out the form below and fax it along with your proof of purchase to **978-762-5223**.

Free Product _____	Cat # _____
YOUR NAME: _____	
TITLE: _____	
INSTITUTION: _____	
DEPT/BLDG: _____	ROOM# _____
ADDRESS 1: _____	
ADDRESS 2: _____	
CITY: _____	
STATE: _____	ZIP _____

To place an order call Fisher Scientific at 1-800-766-7000. For technical support or help in choosing the correct product for your application call Millipore at 1-800-645-5476.



Immobilon is a trademark of Millipore Corporation; Amicon is a registered trademark of Millipore Corporation
 Millipore is a registered trademark of Millipore Corporation. Fisher Scientific is a registered trademark of Fisher Scientific Company LLC

See Reverse For More Offers

FALL FRENZY SPECIAL OFFER



**Buy a Quantum™ & a Q-Gard® Cartridge
Receive a free Millipak® -20 Filter Unit
Now through December 30, 2005**

- **Quantum Cartridge Filters**
- **Q-Gard Purification Packs**
- **MilliPak-20 Express Disposable Filter Units**

Millipore's newest generation of Milli-Q ultrapure water purification systems employs a three-step purification process. Initial purification is provided by a Q-Gard pretreatment pack, which is tailored to your feedwater source. Secondary purification is provided by application-specific Quantum polishing cartridges designed to remove ionic and organic contaminants down to trace levels. Final purification is carried out by a 0.22µm Millipak filter unit for the removal of bacteria and particles above 0.22µm in size.

Ordering Information

QTUM000EX	Quantum EX Cartridge (Organix Resin)	1/Pk
QTUM0001X	Quantum IX Cartridge (Ionex Resin)	1/Pk
QTUM000VX	Quantum VX Cartridge (Volatile Organex Resin)	1/Pk
QGARD00D2	Q-Gard 2 Purification Pak (DI feed water)	1/Pk
QGARD00R1	Q-Gard 1 Purification Pak (Elix RO distilled feed)	1/Pk
MPGP02001	Millipak-20 Filter Units with Express® membrane, 0.22µm, non-sterile	2/Pk

To receive your free Millipak-20 Express Filter Unit, simply fill out the form below and fax it along with your proof of purchase to **978-762-5223**.

NAME: _____
TITLE: _____
INSTITUTION: _____
DEPT/BLDG: _____ ROOM# _____
ADDRESS 1: _____
ADDRESS 2: _____
CITY: _____
STATE: _____ ZIP _____

To place an order call Fisher Scientific at 1-800-766-7000. For technical support or help in choosing the correct product for your application call Millipore at 1-800-645-5476.

